

# MM3

VITAL SIGNS MONITOR



# Mediblu

- Portable Design
- Covidien Quick Temp
- Infrared Ear Thermometer

# MM3 VITAL SIGNS MONITOR

Proven to be a durable and efficient solution for outpatients general wards, the MM3 brings effective vital signs monitoring to most ambulatory and regular healthcares.

- 5.6" high resolution color TFT-LCD screen
- Lightweight & portable design
- User-friendly interface
- Flexible configurations to meet different clinical needs
- Auto/Manual/Continuous/Average BP mode
- Real-time measurements & trend display
- Built-in rechargeable Lithium-ion battery for 8 hours of continuous work
- USB flash disk support
- Barcode scanner support
- Bi-directional with MFM-CMS central monitoring system
- Compatible with PatientCare Viewer PC data management software
- Wired LAN/Wi-Fi support, HL7 support via XML files or MFM-CMS
- Nurse call
- Suitable for adult, pediatric and neonatal patients

**30,000**  
NIBP Measurement  
Review

**72h**  
Trend Review

**800**  
Alarm Review

## Seamless Connectivity



Web Viewer



MFM-CMS



## Mounting Solutions

## Specification

MM3: SPO2+NIBP+COVIDEAN QUICK TEMP



Wall Mount



Rolling Stand



Bed Rail Clamp



**Mediblu**

3511 Silverside Rd, Wilmington, DE, 19810 USA  
Phone: +1 305 406 3931 / [www.mediblumed.com](http://www.mediblumed.com)